


Selection Guide of Hanil Centrifuges

Centrifuge Group		Micro			High Speed, Floor Standing					Multi-pupose, Table-top			Clinical Application, Swing-out			
Model																
		Smart 13	Smart 15 Plus	Smart R17 Plus	Ultra 5.0	Supra R30	Supra R22	Supra R17	Supra R12	Combi R615	Combi 508	Fleta 4	Fleta 5	Component wR	Component R6	Component R12
Max. RPM	Fixed Angle	13,500	15,000	17,000	50,000	30,000	22,000	17,000	12,000	15,000	8,000	4,000	-	-	-	-
	Swing-out	-	-	-	-	-	-	-	-	4,500	4,000	4,000	5,500	4,500	5,000	4,800
Refrigeration		-	-	❄️	❄️	❄️	❄️	❄️	❄️	❄️	-	-	-	❄️	❄️	❄️
Vacuum		-	-	-	⦿	⦿	-	-	-	-	-	-	-	-	-	-



Choices by Sample Container

Container Type	Model	Choice
Micro Plates	Smart 13	
	Smart 15 Plus	
	Smart R17 Plus	
	Ultra 5.0	
	Supra R30	
Micro Tubes	Smart 13	F
	Smart 15 Plus	F
	Smart R17 Plus	F
	Ultra 5.0	F
	Supra R30	F
Tubes	Smart 13	F
	Smart 15 Plus	F
	Smart R17 Plus	F
	Ultra 5.0	F
	Supra R30	F
	Supra R22	F
	Supra R17	F
	Supra R12	F
	Combi R615	F/S
	Combi 508	F/S
	Fleta 4	F/S
	Fleta 5	S
	Component wR	S
	Component R6	S
	Component R12	S
	Syringe or PRP	S
	Blood Bags	S

* F stands for Fixed angle rotor and S for Swing-out rotor.